



The difference
is in the detail.

The difference is in the detail – Acculan 3Ti avoids sterilization not only of the temperature-sensitive battery cells, but also of the equally vulnerable control electronics by integrating it in the battery unit.

A concept which offers you the reliability that you expect and also spares your budget.

Acculan 3Ti dermatome at a glance



■ Acculan 3Ti dermatome GA670

The best dermatome ever from Aesculap. More performance and reliability thanks to the patented motor technology – an instrument you can rely on.

Due to its titanium housing, this dermatome is suitable for use in conjunction with the Eccos storage system for machine reprocessing using alkaline cleaning agents and disinfectants.

- Quick and simple blade changes
- Cutting width max. 78 mm, adjustable by 10 flaps of width 8 mm
- 0.2 – 1.2 mm cutting thickness in 1/10 mm steps, keyless adjustment



■ A clever combination – GA677 charger and GA666 battery

- Intuitive, clear indicator elements
- Automatic battery check
- Charge status display
- Electronically controlled charging to avoid memory effect
- Compatible with predecessor Acculan II through simple-to-replace charging bays
- Innovative battery exchange concept spares the environment and your budget